

ParticleTech Analyzer - Rental Concept

Experience how ParticleTech Analyzer can assist you in your production

- Our rental program combines our **advanced analytical tools** with **expert guidance**
- We tailor **measurement protocol** to your needs, train your team, and provide ongoing support to help you optimize processes and interpret the results
- Our focus is to ensure you gain **actionable insights** during the rental period (minimum of 2 months) and experience first-hand how ParticleTech Analyzer can assist you in your production



ParticleTech Analyzer - Rental Concept

Our rental concept focuses on getting you up and running from day one

Before Rental

- We consult to understand your measurement needs and any potential **operational challenges** or **process optimizations** you aim to use the equipment for during the rental
- We prepare measurement protocols in advance to **ensure efficient** rental use
- You prepare the required workspace (space and power) and relevant samples

Rental Begins

- We install the equipment and load in the pre-made measurement protocols
- We give you a quick introduction to the measurement principles and the software structure
- We provide an **on-site training** of all relevant users*. Each user get **hands-on experience** with the analyzer and the possibility to measure their respective samples

During Rental

- We provide **ongoing technical support** with both operation of the equipment and software
- We provide ongoing **expert guidance** on your measurements to ensure that you obtain the correct data**
- We schedule a follow-up meeting 1-week post-training to review progress and first measurements
- We have **regular check-ins** every 2-3 weeks during the rental for ongoing support

End of Rental

- You pack the equipment back into the designated flight case(s)
- We take care of equipment pickup and return
- We have a follow-up meeting 2 weeks after the end-of-rental to **discuss the outcome** and the results you have obtained

**On-site training takes 1-2 days, depending on the number of users to be trained*

***You retain ownership of all data generated during the rental period*

ParticleTech Analyzer - Rental Concept

Practical preparations before rental

- Bench space (min. 120 cm (L) x 70 cm (D))
- Electrical socket (230 V or 110 V)
- Relevant samples (Recommended quantities per sample: powders 10 g, liquids 50 mL)
- Possibility of sharing files from the rental PC to a OneDrive folder provided by ParticleTech (either via guest Wi-Fi or USB storage)
- Optional: If many users are going to use the instrument, we recommend that each user install the ParticleTech Analyzer software on their respective PC (consult your IT before the rental to inquire whether this is possible)

Personalized support and guidance

- We provide personalized support during the entire rental period via email, phone, WhatsApp and Microsoft Teams
- Furthermore, we also offer analytical guidance to assist you in interpreting the results you obtain during the rental

ParticleTech ApS

Hirsemarken 1

3520 Farum

Denmark

info@particletech.dk

+45 53 54 80 84

